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“Defected soft tissue caused by accidents or trauma”

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"MINT™ which is polydioxanone absorbable surgical suture is intended for correction of defected soft tissue caused by accidents or trauma. MINT™ is inserted under a local anaesthesia in the skin and gently pulled into position. It is the placing of a bared suture under the skin. The suture lift tissues that have defects."

Sikkerhetsprofil beskrives i bruksanvisning som følgende:

Do not re-sterilize. Discard opened, unused MINT™.

To help prevent potential adverse reactions (under Adverse reaction section):

MINT™ only used in hospitals and Physicians should be familiar with surgical procedure a techniques involving absorbable sutures before employing MINT™ for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the in vivo performance (under PERFORMANCES section) when selecting a suture for use in patients.

The safety and effectiveness of MINT™ has not been established for use in fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissue, neural tissue, ophthalmic surgery or for use in microsurgery. Under certain circumstances, notably orthopedic procedures, immobilization by external support may be employed at the discretion of the surgeon.

The use of this suture may be inappropriate in elderly, malnourished or debilitated patients or who are suffering from conditions which may delay wound healing.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, MINT™ may act transiently as a foreign body.

Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention or which may require additional support.

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MINT™ contains bi-directionally oriented barbs to anchor tissues and does not require knots to approximate opposing edges of a wound. Tying knots on the barbed section of the material will damage the barbs and potentially reduce the suture tensile strength and barb

effectiveness. For the bidirectional forces to be created and for the device to function properly, both sides of the MINT™ must be engaged in the tissue. Additionally, when completing placement, an additional backstitch or bite of tissue lateral to the end of the incision is required to lock the device in place. ·Avoid contacting MINT™ with other equipment (e.g. surgical gauzes, drapes etc) in the surgical field to prevent ensnaring of the barbs. If the barbs catch, carefully pull the material in the opposite direction to disengage it from the barbs. ·Care should be taken to avoid damage when handling. Avoid crushing or crimping the suture by sterilized surgical instruments, such as sterilized needle holders and forceps. Do not pull the MINT™ out of the package by the needles as this can cause the barbs to catch on one another. Do not attempt to remove memory in the polymer by running fingers down the suture material as this can damage the barbs. ·Infections, erythema, foreign body reactions, transient inflammatory reactions and in rare instances wound dehiscence are typical or foreseeable risks associated with any suture, and hence are also potential complications associated with MINT™. User must use with sterilized medical needles that consists of biocompatibility materials (medical grade stainless steel). In order to minimize the above risks. ·When using MINT™ subcutaneously, the device should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption.

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MINT™ is not to be used where prolonged approximation of tissues under stress is required and is not be used in conjunction with or for fixation of prosthetic devices (e.g. heart valves or synthetic grafts) that are non-absorbable in nature.

Produsent: HansBiomed Corp. 64, Yuseong-daero 1628beon-gil, Yuseong-gu, Daejeon, 34054, Republic of Korea

Autorisert representant: EMERGO EUROPE, Prinsessegracht 20
2514 AP The Hague, The Netherlands

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“Results of implantation studies in rats indicate that approximately 88% of the original tensile strength remains at 2 weeks post implantation, approximately 70% remains at 4 weeks, and approximately 45% remains at 6 weeks. Data from implantation studies show that absorption

is minimal until about 120 days after implantation, that significant suture absorption occurs between 120-180 days after implantation, and that absorption is essentially complete”.

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«This revolutionary technology allows instant yet profound results while minimizing pain & trauma - allowing patients to benefit from a shorter recovery time and providing an alternative for those who may not be ready for surgery or request a non-surgical option.» (Publisert 13.12.2021, vedlegg 1 side 1)

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